

Neonovo IV Cannula						
Ref. Code	Gauge	Catheter Length (mm)	Catheter OD (mm)	Flow Rate (ml/min)	Catheter Material	Packaging
10615	26 ●	19	0.60	17	FEP	100 each/Duplex & 10 Duplex/Case
10616	26 ●	19	0.60	17	FEP	100 each/Duplex & 10 Duplex/Case
10617	26 ●	19	0.60	17	FEP	100 each/Duplex & 10 Duplex/Case
10614	26 ●	19	0.60	17	FEP	100 each/Duplex & 10 Duplex/Case
10619	26 ●	19	0.60	17	FEP	100 each/Duplex & 10 Duplex/Case
10108	26 ●	19	0.60	17	FEP	100 each/Duplex & 10 Duplex/Case

Packaging Dimension	
Unit Pack Size	96mm X 27mm
Duplex Box Size	170mm X 100mm X 150mm
Case Size	525mm X 320mm X 188mm