

Nuovo Safety Set IV Cannula											
Gauge	Catheter Length (mm)	Catheter OD (mm)	Flow Rate (ml/min)	Catheter Material	Product Code (Y-site)		Double Luer Lock	Straight Luer Lock	Stop Cock	Stop Cock with NFC*	Packaging
					Single NFC**		Double NFC*				
16 	32	1.7	125	PUR	11988	11980	11922	11930	11945	11936	20 each/Duplex & 10 Duplex/Case
18 	45	1.3	79	PUR	11990	11982	11954	11955	11946	11937	20 each/Duplex & 10 Duplex/Case
18 	32	1.3	84	PUR	11989	11981	11923	11931	11947	11938	20 each/Duplex & 10 Duplex/Case
20 	32	1.1	58	PUR	11991	11983	11924	11932	11948	11939	20 each/Duplex & 10 Duplex/Case
20 	25	1.1	61	PUR	11786	11785	11782	11783	11784	11788	20 each/Duplex & 10 Duplex/Case
22 	25	0.9	33	PUR	11992	11984	11925	11933	11949	11940	20 each/Duplex & 10 Duplex/Case
24 	19	0.7	18	PUR	11993	11985	11926	11934	11950	11941	20 each/Duplex & 10 Duplex/Case
26 	19	0.6	15	PUR	11994	11986	11927	11935	11951	11942	20 each/Duplex & 10 Duplex/Case

Packaging Dimension	
Unit Pack Size	128mm X 78mm
Duplex Box Size	160mm X 132mm X 132mm
Case Size	685mm X 275mm X 180mm