

Polycan IV Cannula						
Ref. Code	Gauge	Catheter Length (mm)	Catheter OD (mm)	Flow Rate (ml/min)	Catheter Material	Packaging
10100	14 	45	2.10	305	FEP	100 each/Duplex & 10 Duplex/Case
10101	16 	45	1.70	200	FEP	100 each/Duplex & 10 Duplex/Case
10102	17 	45	1.50	142	FEP	100 each/Duplex & 10 Duplex/Case
10103	18 	45	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
10104	18 	32	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
10105	20 	32	1.10	61	FEP	100 each/Duplex & 10 Duplex/Case
10106	22 	25	0.90	36	FEP	100 each/Duplex & 10 Duplex/Case
10107	24 	19	0.70	23	FEP	100 each/Duplex & 10 Duplex/Case
10108	26 	19	0.60	17	FEP	100 each/Duplex & 10 Duplex/Case
10136	14 	45	2.10	305	PUR	100 each/Duplex & 10 Duplex/Case
10137	16 	45	1.70	200	PUR	100 each/Duplex & 10 Duplex/Case
10138	17 	45	1.50	142	PUR	100 each/Duplex & 10 Duplex/Case
10139	18 	45	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
10140	18 	32	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
10141	20 	32	1.10	61	PUR	100 each/Duplex & 10 Duplex/Case
10142	22 	25	0.90	36	PUR	100 each/Duplex & 10 Duplex/Case
10143	24 	19	0.70	23	PUR	100 each/Duplex & 10 Duplex/Case
10144	26 	19	0.60	17	PUR	100 each/Duplex & 10 Duplex/Case

Packaging Dimension	
Unit Pack Size	114mm X 27mm
Duplex Box Size	170mm X 117mm X 150mm
Case Size	605mm X 320mm X 188mm