

Polyflon IV Cannula						
Ref. Code	Gauge	Catheter Length (mm)	Catheter OD (mm)	Flow Rate (ml/min)	Catheter Material	Packaging
10001	14 	45	2.10	305	FEP	100 each/Duplex & 10 Duplex/Case
10002	16 	45	1.70	200	FEP	100 each/Duplex & 10 Duplex/Case
10003	17 	45	1.50	142	FEP	100 each/Duplex & 10 Duplex/Case
10004	18 	45	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
10005	18 	32	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
10006	20 	32	1.10	61	FEP	100 each/Duplex & 10 Duplex/Case
10007	22 	25	0.90	36	FEP	100 each/Duplex & 10 Duplex/Case
10008	24 	19	0.70	23	FEP	100 each/Duplex & 10 Duplex/Case
10009	26 	19	0.60	17	FEP	100 each/Duplex & 10 Duplex/Case
10055	14 	45	2.10	305	PUR	100 each/Duplex & 10 Duplex/Case
10056	16 	45	1.70	200	PUR	100 each/Duplex & 10 Duplex/Case
10057	17 	45	1.50	142	PUR	100 each/Duplex & 10 Duplex/Case
10058	18 	45	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
10059	18 	32	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
10060	20 	32	1.10	61	PUR	100 each/Duplex & 10 Duplex/Case
10061	22 	25	0.90	36	PUR	100 each/Duplex & 10 Duplex/Case
10062	24 	19	0.70	23	PUR	100 each/Duplex & 10 Duplex/Case
10063	26 	19	0.60	17	PUR	100 each/Duplex & 10 Duplex/Case

Packaging Dimension	
Unit Pack Size	114mm X 27mm
Duplex Box Size	170mm X 117mm X 150mm
Case Size	605mm X 320mm X 188mm