

Polyneo IV Cannula						
Ref. Code	Gauge	Catheter Length (mm)	Catheter OD (mm)	Flow Rate (ml/min)	Catheter Material	Packaging
10600	24 ●	19	0.70	23	FEP	100 each/Duplex & 10 Duplex/Case
10601	24 ●	19	0.70	23	FEP	100 each/Duplex & 10 Duplex/Case
10602	24 ●	19	0.70	23	FEP	100 each/Duplex & 10 Duplex/Case
10603	24 ●	19	0.70	23	FEP	100 each/Duplex & 10 Duplex/Case
10604	24 ●	19	0.70	23	FEP	100 each/Duplex & 10 Duplex/Case
10612	24 ●	19	0.70	23	FEP	100 each/Duplex & 10 Duplex/Case
10610	24 ●	19	0.70	23	PUR	100 each/Duplex & 10 Duplex/Case
10611	24 ●	19	0.70	23	PUR	100 each/Duplex & 10 Duplex/Case
10613	24 ●	19	0.70	23	PUR	100 each/Duplex & 10 Duplex/Case

Packaging Dimension	
Unit Pack Size	96mm X 27mm
Duplex Box Size	170mm X 100mm X 150mm
Case Size	525mm X 320mm X 188mm