

Polysafety IV Cannula						
Ref. Code	Gauge	Catheter Length (mm)	Catheter OD (mm)	Flow Rate (ml/min)	Catheter Material	Packaging
11001	14 	45	2.10	305	FEP	100 each/Duplex & 10 Duplex/Case
11002	16 	45	1.70	200	FEP	100 each/Duplex & 10 Duplex/Case
11003	17 	45	1.50	142	FEP	100 each/Duplex & 10 Duplex/Case
11004	18 	45	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
11005	18 	32	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
11006	20 	32	1.10	61	FEP	100 each/Duplex & 10 Duplex/Case
11007	22 	25	0.90	36	FEP	100 each/Duplex & 10 Duplex/Case
11008	24 	19	0.70	23	FEP	100 each/Duplex & 10 Duplex/Case
11009	26 	19	0.60	17	FEP	100 each/Duplex & 10 Duplex/Case
11010	14 	45	2.10	305	PUR	100 each/Duplex & 10 Duplex/Case
11011	16 	45	1.70	200	PUR	100 each/Duplex & 10 Duplex/Case
11012	17 	45	1.50	142	PUR	100 each/Duplex & 10 Duplex/Case
11013	18 	45	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
11014	18 	32	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
11015	20 	32	1.10	61	PUR	100 each/Duplex & 10 Duplex/Case
11016	22 	25	0.90	36	PUR	100 each/Duplex & 10 Duplex/Case
11017	24 	19	0.70	23	PUR	100 each/Duplex & 10 Duplex/Case
11018	26 	19	0.60	17	PUR	100 each/Duplex & 10 Duplex/Case

Packaging Dimension	
Unit Pack Size	138mm X 32mm
Duplex Box Size	170mm X 135mm X 105mm
Case Size	540mm X 300mm X 185mm