

Polyspike		
Product Code	Product Description	Units / Duplex Box
13068	Polyspike -Air vented spike with green snap fit cap, soft blister	100
13490	Polyspike- Air vented spike with Blue snap fit cap, fluid filter, soft blister	100
13491	Polyspike- Air vented chemo spike with Red snap fit cap, fluid filter, soft blister	100
13492	Polyspike -Air vented spike with needle free valve, soft blister	100
13493	Polyspike- Air vented spike with needle free valve, fluid filer, soft blister	100
13494	Polyspike- Air vented chemo spike with needle free valve, fluid filter, soft blister	100
13495	Polyspike -Air vented spike with green snap fit cap (air filter 0.1 µm), soft blister	100
13496	Polyspike- Air vented spike with Blue snap fit cap (air filter 0.1 µm), fluid filter (fluid filter 5 µm), soft blister	100