

Polywin IV Cannula						
Ref. Code	Gauge	Catheter Length (mm)	Catheter OD (mm)	Flow Rate (ml/min)	Catheter Material	Packaging
10210	14 	45	2.10	305	FEP	100 each/Duplex & 10 Duplex/Case
10211	16 	45	1.70	200	FEP	100 each/Duplex & 10 Duplex/Case
10212	17 	45	1.50	142	FEP	100 each/Duplex & 10 Duplex/Case
10213	18 	45	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
10214	18 	32	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
10215	20 	32	1.10	61	FEP	100 each/Duplex & 10 Duplex/Case
10216	22 	25	0.90	36	FEP	100 each/Duplex & 10 Duplex/Case
10217	24 	19	0.70	23	FEP	100 each/Duplex & 10 Duplex/Case
10218	26 	19	0.60	17	FEP	100 each/Duplex & 10 Duplex/Case
10250	14 	45	2.10	305	PUR	100 each/Duplex & 10 Duplex/Case
10251	16 	45	1.70	200	PUR	100 each/Duplex & 10 Duplex/Case
10252	17 	45	1.50	142	PUR	100 each/Duplex & 10 Duplex/Case
10253	18 	45	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
10254	18 	32	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
10255	20 	32	1.10	61	PUR	100 each/Duplex & 10 Duplex/Case
10256	22 	25	0.90	36	PUR	100 each/Duplex & 10 Duplex/Case
10257	24 	19	0.70	23	PUR	100 each/Duplex & 10 Duplex/Case
10258	26 	19	0.60	17	PUR	100 each/Duplex & 10 Duplex/Case

Packaging Dimension	
Unit Pack Size	134mm X 27mm
Duplex Box Size	170mm X 135mm X 105mm
Case Size	540mm X 300mm X 185mm