

Polywin Plus IV Cannula

Ref. Code	Gauge	Catheter Length (mm)	Catheter OD (mm)	Flow Rate (ml/min)	Catheter Material	Packaging
10450	14 	45	2.10	305	PUR	100 each/Duplex & 10 Duplex/Case
10451	16 	45	1.70	200	PUR	100 each/Duplex & 10 Duplex/Case
10452	17 	45	1.50	142	PUR	100 each/Duplex & 10 Duplex/Case
10453	18 	45	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
10454	18 	32	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
10455	20 	32	1.10	61	PUR	100 each/Duplex & 10 Duplex/Case
10456	22 	25	0.90	36	PUR	100 each/Duplex & 10 Duplex/Case
10457	24 	19	0.70	23	PUR	100 each/Duplex & 10 Duplex/Case
10458	26 	19	0.60	17	PUR	100 each/Duplex & 10 Duplex/Case
10460	14 	45	2.10	305	FEP	100 each/Duplex & 10 Duplex/Case
10461	16 	45	1.70	200	FEP	100 each/Duplex & 10 Duplex/Case
10462	17 	45	1.50	142	FEP	100 each/Duplex & 10 Duplex/Case
10463	18 	45	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
10464	18 	32	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
10465	20 	32	1.10	61	FEP	100 each/Duplex & 10 Duplex/Case
10466	22 	25	0.90	36	FEP	100 each/Duplex & 10 Duplex/Case
10467	24 	19	0.70	23	FEP	100 each/Duplex & 10 Duplex/Case
10468	26 	19	0.60	17	FEP	100 each/Duplex & 10 Duplex/Case

Packaging Dimension

Unit Pack Size	155mm X 32mm
Duplex Box Size	170mm X 135mm X 105mm
Case Size	540mm X 300mm X 185mm