

Polywin Safety IV Cannula

Ref. Code	Gauge	Catheter Length (mm)	Catheter OD (mm)	Flow Rate (ml/min)	Catheter Material	Packaging
11090	14 	45	2.10	305	FEP	100 each/Duplex & 10 Duplex/Case
11091	16 	45	1.70	200	FEP	100 each/Duplex & 10 Duplex/Case
11092	17 	45	1.50	142	FEP	100 each/Duplex & 10 Duplex/Case
11093	18 	45	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
11094	18 	32	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
11095	20 	32	1.10	61	FEP	100 each/Duplex & 10 Duplex/Case
11096	22 	25	0.90	36	FEP	100 each/Duplex & 10 Duplex/Case
11097	24 	19	0.70	23	FEP	100 each/Duplex & 10 Duplex/Case
11098	26 	19	0.60	17	FEP	100 each/Duplex & 10 Duplex/Case
11100	14 	45	2.10	305	PUR	100 each/Duplex & 10 Duplex/Case
11101	16 	45	1.70	200	PUR	100 each/Duplex & 10 Duplex/Case
11102	17 	45	1.50	142	PUR	100 each/Duplex & 10 Duplex/Case
11103	18 	45	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
11104	18 	32	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
11105	20 	32	1.10	61	PUR	100 each/Duplex & 10 Duplex/Case
11106	22 	25	0.90	36	PUR	100 each/Duplex & 10 Duplex/Case
11107	24 	19	0.70	23	PUR	100 each/Duplex & 10 Duplex/Case
11108	26 	19	0.60	17	PUR	100 each/Duplex & 10 Duplex/Case

Packaging Dimension

Unit Pack Size	134mm X 27mm
Duplex Box Size	170mm X 135mm X 105mm
Case Size	540mm X 300mm X 185mm