

Haemocheck (Coagulation (Citrate) Tubes (Single wall)						
---	--	--	--	--	--	--

Product Code	Draw Volume	Cap Color		Tube Size (In mm)	Additive	Units / Duplex Box
81020	1.8 ml	Light Blue	●	13 x 75	3.2% Buffered Sodium Citrate 9NC	100
81030	2.7 ml	Light Blue	●	13 x 75	3.2% Buffered Sodium Citrate 9NC	100
81040	4.5 ml	Light Blue	●	13 x 75	3.2% Buffered Sodium Citrate 9NC	100

Packaging Dimension	
Unit Pack Size	114mm X 27mm
Duplex Box Size	NA
Case Size	550mm X 465mm X 195mm