

Haemocheck (Lithium Heparin Tubes)						
Product Code	Draw Volume	Cap Color		Tube Size (In mm)	Additive	Units / Duplex Box
88030	3 ml	Green	●	13 x 75	Lithium Heparin	100
88040	4 ml	Green	●	13 x 75	Lithium Heparin	100
88050	5 ml	Green	●	13 x 100	Lithium Heparin	100
88060	6 ml	Green	●	13 x 100	Lithium Heparin	100
88070	7 ml	Green	●	13 x 100	Lithium Heparin	100
88100	10 ml	Green	●	16 x 100	Lithium Heparin	100

Packaging Dimension	
Unit Pack Size	114mm X 27mm
Duplex Box Size	NA
Case Size	550mm X 465mm X 195mm