

Polycan Safety IV Cannula

Ref. Code	Gauge	Catheter Length (mm)	Catheter OD (mm)	Flow Rate (ml/min)	Catheter Material	Packaging
11050	14 	45	2.10	305	FEP	100 each/Duplex & 10 Duplex/Case
11051	16 	45	1.70	200	FEP	100 each/Duplex & 10 Duplex/Case
11052	17 	45	1.50	142	FEP	100 each/Duplex & 10 Duplex/Case
11053	18 	45	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
11054	18 	32	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
11055	20 	32	1.10	61	FEP	100 each/Duplex & 10 Duplex/Case
11056	22 	25	0.90	36	FEP	100 each/Duplex & 10 Duplex/Case
11057	24 	19	0.70	23	FEP	100 each/Duplex & 10 Duplex/Case
11058	26 	19	0.60	17	FEP	100 each/Duplex & 10 Duplex/Case
11060	14 	45	2.10	305	PUR	100 each/Duplex & 10 Duplex/Case
11061	16 	45	1.70	200	PUR	100 each/Duplex & 10 Duplex/Case
11062	17 	45	1.50	142	PUR	100 each/Duplex & 10 Duplex/Case
11063	18 	45	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
11064	18 	32	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
11065	20 	32	1.10	61	PUR	100 each/Duplex & 10 Duplex/Case
11066	22 	25	0.90	36	PUR	100 each/Duplex & 10 Duplex/Case
11067	24 	19	0.70	23	PUR	100 each/Duplex & 10 Duplex/Case
11068	26 	19	0.60	17	PUR	100 each/Duplex & 10 Duplex/Case

Packaging Dimension

Unit Pack Size	134mm X 27mm
Duplex Box Size	170mm X 135mm X 105mm
Case Size	540mm X 300mm X 185mm