

Polywin Plus Safety IV Cannula

Ref. Code	Gauge	Catheter Length (mm)	Catheter OD (mm)	Flow Rate (ml/min)	Catheter Material	Packaging
11400	14 	45	2.10	305	FEP	100 each/Duplex & 10 Duplex/Case
11401	16 	45	1.70	200	FEP	100 each/Duplex & 10 Duplex/Case
11402	17 	45	1.50	142	FEP	100 each/Duplex & 10 Duplex/Case
11403	18 	45	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
11404	18 	32	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
11405	20 	32	1.10	61	FEP	100 each/Duplex & 10 Duplex/Case
11406	22 	25	0.90	36	FEP	100 each/Duplex & 10 Duplex/Case
11407	24 	19	0.70	23	FEP	100 each/Duplex & 10 Duplex/Case
11408	26 	19	0.60	17	FEP	100 each/Duplex & 10 Duplex/Case
11420	14 	45	2.10	305	PUR	100 each/Duplex & 10 Duplex/Case
11421	16 	45	1.70	200	PUR	100 each/Duplex & 10 Duplex/Case
11422	17 	45	1.50	142	PUR	100 each/Duplex & 10 Duplex/Case
11423	18 	45	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
11424	18 	32	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
11425	20 	32	1.10	61	PUR	100 each/Duplex & 10 Duplex/Case
11426	22 	25	0.90	36	PUR	100 each/Duplex & 10 Duplex/Case
11427	24 	19	0.70	23	PUR	100 each/Duplex & 10 Duplex/Case
11428	26 	19	0.60	17	PUR	100 each/Duplex & 10 Duplex/Case

Packaging Dimension

Unit Pack Size	134mm X 27mm
Duplex Box Size	170mm X 135mm X 105mm
Case Size	540mm X 300mm X 185mm